



CHECKING OCCLUSION USING THE 2-STEP METHOD

To visualize contacts on saliva moistened surfaces, we recommend using a progression of test materials from thicker papers to thinner foils.



200 μ

100 μ

12 μ



Articulating Papers

Arti-Fol® Foil

Our *PAPERS* absorb saliva and create the perfect foundation for the pinpoint accurate markings of our *FOILS*!

ALL OF OUR PRODUCTS ARE AVAILABLE IN ADDITIONAL SHAPES, THICKNESSES, AND/OR COLORS!

Bausch
WE MAKE OCCLUSION VISIBLE