


1. Protocol: Occlusal Contacts

© Dr. Jean Bausch GmbH & Co. KG 2022


To mark position of initial occlusal contacts before adjustment.

Date:

Quadrant 1



Quadrant 2




Doc.No.:


Patient:

Dentist:

Quadrant 4



Quadrant 3



Notes:

2. Protocol: Shimstock Hold

To archive position of where Shimstock is held or pulled through

Put your patient in an upright position. The teeth should be in maximum intercuspitation. The patient bites on it only slightly, with low muscle tone. Check tooth for tooth, and enter the result into the sheet.

18	17	16	15	14	13	12	11	21	22	23	24	25	26	27	28
48	47	46	45	44	43	42	41	31	32	33	34	35	36	37	38

X = hold ; 0 = pulled through (no tactile contact)

Recommended Test Materials:

Bausch Articulating Papers ø 100 micron Ref-No.: BK 51, BK 52, BK 53, BK 54....

Bausch Articulating Papers ø 40 micron Ref-No.: BK 09, BK 10, BK 80, BK 81....

Bausch Arti-Fol Metallic (Shimstock) ø 12 micron Ref-No.: BK 30, 31, 32, 33

Bausch Arti-Fol Metallic (Shimstock) precut strips ø 12 micron Ref-No.: BK 35, BK 38

Bausch Arti-Fol Forceps BK 132 + Arti-Grip Silicon sleeves BK 146

Bausch
WE MAKE OCCLUSION VISIBLE

Dr. Jean Bausch GmbH & Co. KG | Oskar-Schindler-Str. 4 | D-50769 Köln - Germany
Phone: +49-221-70936-0 | Fax: ++49-221-70936-66
info@bauschdental.de | www.bauschdental.de