

ARTI-CHECK 40 μ ARTICULATING PAPER

Bausch Articulating Papers provide the highest quality and most accurate results in every procedure. *How do we accomplish this?* Our marking material is embedded in special microfleece paper through a patented process resulting in a progressive color transfer that helps you achieve accurate occlusion—every time.

THE PRODUCT LINE



BK13 & BK14
16mm Wide



BK15 & BK16
22mm Wide



BK09, BK10, & BK80
200 Strips



BK61, BK62, & BK63
200 Strips


BK11 & BK12
100 Sheets




WHAT MAKES OUR PRODUCT BETTER?

- Works well on gold, ceramic, polished metal, or acrylic
- Excels in both wet and dry applications
- Color on *both* sides
- Precise marking to avoid false marks
- Different versions can be used with or without forceps

Our papers are extremely thin and tear resistant while still providing highly visible markings.



Even the slightest masticatory pressure can cause the capsules to burst and thus release the distinctly visible color.



- Horseshoe-shaped Articulating Papers are especially useful for patients who tend to bite unilaterally during the occlusion test. The dentist can immediately detect the preferred side of the mouth.

**BK17, BK18,
& BK81**
150 Sheets



Bausch
WE MAKE OCCLUSION VISIBLE

More info at
BAUSCHPAPER.COM